

SD_Study Data Files Scanner and Analyzer

In the drug registration process, SD files (short for Study Data files) refer to electronic data submissions that contain structured datasets from nonclinical and clinical studies. These datasets are prepared following standardized data models to allow regulatory agencies (like the FDA, EMA, or local authorities) to efficiently review, validate, and analyze the study results.

Definition

SD files = Study Data Files — electronic data packages containing the raw and tabulated data from: **Nonclinical studies** (toxicology, pharmacology) **Clinical studies** (efficacy, safety, pharmacokinetics, etc.)

CDISC standards

- SDTM (Study Data Tabulation Model) – standardized structure for tabulated data.
- ADAM (Analysis Data Model) – datasets used for statistical analysis.
- SEND (Standard for Exchange of Nonclinical Data) – for preclinical (animal) data.

This notebook scans the current directory for files and analyzes specific file types using appropriate libraries:

- **.sdf files:** Structure Data File (Chemoinformatics) Plain text molecular data format used to store 3D structures, bonds, properties, and metadata - RDKit for chemical structure analysis
- **.xpt files:** SAS XPORT (SDTM Clinical Datasets) Contains clinical study data such as Demographics (DM), Adverse Events (AE), and Laboratory Results (LB) - pyreadstat for SAS transport files
- **.asn1 files:** Textual representation of Abstract Syntax Notation One, used in regulatory metadata, pharma labeling, and bioinformatics standards - Biopython for ASN.1 files

In [1]: # Install required packages if not already installed
Uncomment the following lines if you need to install packages

```
!pip install -q rdkit-pypi  
!pip install -q pandas  
!pip install -q biopython
```

In [2]:
import os
import glob
from pathlib import Path

```
# Scan current directory for files  
current_dir = Path('.')  
all_files = list(current_dir.glob('*'))  
files_only = [f for f in all_files if f.is_file()]
```

```
print(f"Found {len(files_only)} files in the current directory and subdirectories")
```

```
# Group files by extension
```

```
file_extensions = {}  
for file_path in files_only:  
    ext = file_path.suffix.lower()  
    if ext not in file_extensions:  
        file_extensions[ext] = []  
    file_extensions[ext].append(file_path)
```

```
print(f"Files by extension")
```

```
for ext, files in file_extensions.items():  
    print(f"  {ext} {len(files)} files")
```

```
    for file in files[:5]: # Show first 5 files per extension
```

```
        print(f"    {file}")
```

```
    if len(files) > 5:  
        print(f"    ... and {len(files) - 5} more")
```

```
Found 25 files in the current directory and subdirectories
```

```
Files by extension:
```

```
5 files
```

```
- DS_Store
```

```
- m53-DS_Store
```

```
- m53-clin-stud-reports/.DS_Store
```

```
- m53-clin-stud-reports/study1234/.DS_Store
```

```
- m53-clin-stud-reports/study1234/datasets/.DS_Store
```

```
- SD_Study-Data-files.pywb
```

```
- xlsx 1 files
```

```
- sample-data-summary.xlsx
```

```
json: 1 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/COMPOUND_CID_197365.json
```

```
xml: 3 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/strf.xml
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/define.xml
```

```
xpt: 3 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/lb.xpt
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/ae.xpt
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/dm.xpt
```

```
pdf: 1 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/Conformer3D_COMPOUND_CID_197366.pdf
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/Structure2D_COMPOUND_CID_197365.pdf
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/assessment.pdf
```

```
pdf: 1 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/annotated-crf.pdf
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/study1234-clin-report.pdf
```

```
- xsl: 3 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/compound.xsl
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/Structure2D_COMPOUND_CID_197365.xsl
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/dm.xsl
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/ae.xsl
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/lb.xsl
```

```
csv: 3 files
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/dm.csv
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/ae.csv
```

```
- m53-clin-stud-reports/study1234/datasets/datasets/lb.csv
```

SDF Files Analysis (RDKit)

In [3]: # Analyze SDF files using RDKit

```
try:  
    from rdkit import Chem  
    from rdkit.Chem import Descriptors
```

```
    sdf_files = file_extensions.get('.sdf', [])
```

```
    if sdf_files:  
        print(f"Found {len(sdf_files)} SDF files")
```

```
        for sdf_file in sdf_files:  
            print(f"Analyzing {sdf_file}:")
```

```
            # Read SDF file
```

```
            suppl = Chem.SDMolSupplier(str(sdf_file))
```

```
            molecules = [mol for mol in suppl if mol is not None]
```

```
            print(f"  - Number of molecules: {len(molecules)}")
```

```
            if molecules:  
                # Analyze first molecule as example
```

```
                mol = molecules[0]
```

```
                print(f"  - First molecule: {mol.GetProp('_Name') if mol.HasProp('_Name') else 'Unnamed'}")
```

```
                print(f"  - Molecular weight: {Descriptors.MolWt(mol):.2f}")
```

```
                print(f"  - Number of atoms: {mol.GetNumAtoms()}")
```

```
                print(f"  - Number of bonds: {mol.GetNumBonds()}")
```

```
                print(f"  - SMILES: {Chem.MolToSmiles(mol)}")
```

```
            else:  
                print("No SDF files found")
```

```
except ImportError:  
    print("RDKit not installed. Install with: pip install rdkit-pypi")
```

```
except Exception as e:  
    print(f"Error analyzing SDF files: {e}")
```

```
Found 3 SDF files
```

```
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/compound.sdf:
```

```
- Number of molecules: 0
```

```
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/Structure2D_COMPOUND_CID_197365.sdf:
```

```
- Number of molecules: 1
```

```
- First molecule: 197365
```

```
- Molecular weight: 699.99
```

```
- Number of atoms: 46
```

```
- Number of bonds: 46
```

```
- SMILES: O=C(CN2C=CC=NC(=C2C)C(C)C)C(C)C
```

```
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/Conformer3D_COMPOUND_CID_197366.sdf:
```

```
- Number of molecules: 1
```

```
- First molecule: 197366
```

```
- Molecular weight: 598.68
```

```
- Number of atoms: 40
```

```
- Number of bonds: 46
```

```
- SMILES: O=C(CN2C=CC=NC(=C2C)C(C)C)C(C)C
```

```
16:12:02[1]: ERROR: Atom line too short: '-0.0015 1.2095 0.0000 C' on line 5
```

```
16:12:02[1]: ERROR: moving to the beginning of the next molecule
```

XPT Files Analysis (pyreadstat)

In [4]: # Analyze XPT files using pyreadstat

```
try:  
    import pyreadstat
```

```
    xpt_files = file_extensions.get('.xpt', [])
```

```
    if xpt_files:  
        print(f"Found {len(xpt_files)} XPT files")
```

```
        for xpt_file in xpt_files:  
            print(f"Analyzing {xpt_file}:")
```

```
            try:  
                # Read XPT file using pyreadstat
```

```
                df, meta = pyreadstat.read_xport(str(xpt_file))
```

```
                print(f"  - Shape: {df.shape}")
```

```
                print(f"  - Columns: {list(df.columns)}")
```

```
                print(f"  - First 5 rows:{df.head()}\")
```

```
            except Exception as file_error:  
                print(f"  - Error reading {xpt_file}: {file_error}")
```

```
        else:  
            print("No XPT files found")
```

```
except ImportError as e:  
    print(f"Required packages not installed: {e}")
```

```
    print("Install with pip install pyreadstat")
```

```
except Exception as e:  
    print(f"Error analyzing XPT files: {e}")
```

```
...
```

```
Found 3 XPT files
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/lb.xpt:
- Shape: (2, 5)
- Columns: ['STUDYID', 'SUBJID', 'LBTEST', 'LBSTRESN', 'LBSTRESU']
- First 5 rows:
STUDYID SUBJID LBTEST LBSTRESN LBSTRESU
0 ABC123 SUBJ001 Hemoglobin 14.2 g/dl
1 ABC123 SUBJ002 Glucose 89.8 mg/dl
- Summary statistics:
LBSTRESN
count 2.000000
mean 51.100000
std 52.18448
min 14.200000
25% 32.650000
50% 51.100000
75% 69.550000
max 88.800000
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/ae.xpt:
- Shape: (2, 5)
- Columns: ['STUDYID', 'SUBJID', 'AETERM', 'AESEV', 'AEREL']
- First 5 rows:
STUDYID SUBJID AETERM AESEV AEREL
0 ABC123 SUBJ001 Headache MILD RELATED
1 ABC123 SUBJ002 Nausea MODERATE UNRELATED
- Summary statistics:
AETERM AESEV AEREL
count 2 2 2
unique 1 2 2
top ABC123 SUBJ001 Headache MILD RELATED
freq 1 1 1
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/dm.xpt:
- Shape: (2, 4)
- Columns: ['STUDYID', 'SUBJID', 'AGE', 'SEX']
- First 5 rows:
STUDYID SUBJID AGE SEX
0 ABC123 SUBJ001 45.0 M
1 ABC123 SUBJ002 52.0 F
- Summary statistics:
AGE
count 2.000000
mean 48.500000
std 4.949747
min 45.000000
25% 46.750000
50% 48.500000
75% 50.250000
max 52.000000
```

ASNT Files Analysis (Biopython)

```
In [5]: # Analyze ASNT files using Biopython and XML parsing
try:
    from Bio import SeqIO
    import xml.etree.ElementTree as ET

    def read_asnt(asnt_file):
        with open(asnt_file, 'r') as f:
            first_line = f.readline()
            if first_line.startswith('<?xml'):
                if f.readline().startswith('<!DOCTYPE'):
                    tree = ET.parse(asnt_file)
                    root = tree.getroot()
                    return [root] # Return as list with one element
                else:
                    # Try as ASNT sequence
                    records = list(SeqIO.parse(asnt_file, "genbank"))
                    return records
            asnt_files = file_extensions.get('.asnt', [])
            if asnt_files:
                print(f"Found {len(asnt_files)} ASNT files")
            for asnt_file in asnt_files:
                print(f"Analyzing {asnt_file};")

            try:
                records = read_asnt(asnt_file)
                print(f" - Number of records: {len(records)}")

                for i, record in enumerate(records[3:]): # Show first 3 records
                    if hasattr(record, 'id'):
                        print(f" - Record {i+1}: {record.id}")
                        print(f"   - Record {i+1}: {record.id}")
                        print(f"   - Description: {record.description}")
                        print(f"   - Sequence length: {len(record.seq)}")
                        print(f"   - Sequence type: {record.seq.alphabet}")
                    else:
                        # XML element
                        print(f" - XML Root: {record.tag}")
                        for child in record:
                            print(f"   - {child.tag}: {child.text}")
                print("...")

            except Exception as e:
                print(f" - Error reading {asnt_file}: {e}")
                # Fallback: read as text
                try:
                    with open(asnt_file, 'r') as f:
                        content = f.read()
                        print(f" - File size: {len(content)} characters")
                        print(f" - First 500 characters:\n{content[:500]}...")
                except Exception as e2:
                    print(f" - Error reading as text: {e2}")
            else:
                print("No ASNT files found")
        except IOError:
            print("Biopython not installed. Install with: pip install biopython")
        except Exception as e:
            print(f"Error analyzing ASNT files: {e}")
Found 3 ASNT files

Analyzing m53-clin-stud-reports/study1234/datasets/datasets/Conformer3D_COMPOUND_CID_197366.asnt:
- Number of records: 0

Analyzing m53-clin-stud-reports/study1234/datasets/datasets/Structure2D_COMPOUND_CID_197365.asnt:
- Number of records: 0
--
```

```
Analyzing m53-clin-stud-reports/study1234/datasets/datasets/assessment.asnt:
- Number of records: 1
- XML Root: AssessmentTemplate
- Reviewer: Dr. Smith
- AssessmentDate: 2025-11-05
Findings:
-- Recommendation: Approved for next review phase
```

Summary

```
In [6]: # Summary of analysis
# Summary of analysis
# Summary of analysis
print("SUMMARY ===")
print(f"Total files scanned: {len(files_only)}")
print(f"File extensions found: {len(file_extensions)}")
print(f"File extensions found: {len(file_extensions)}")

# Check specific file types
special_types = ['.sdf', '.xpt', '.asnt']
for ext in special_types:
    count = len(file_extensions.get(ext, []))
    print(f"Ext {ext}: {count}")

print("\nAnalysis complete!")

# THIS SUMMARY ==
Total files scanned: 20
File extensions found: 10
.SDF files: 3
.XPT files: 3
.ASNT files: 3
Analysis complete!
```

Enhanced Comprehensive File Summary Table with Detailed Content Columns

```
In [7]: # Create enhanced summary table with detailed content columns
import pandas as pd

# Helper functions to extract detailed file information
def analyze_sdf(file_path):
    """Analyze SDF file and return detailed information."""
    try:
        from rdkit import Chem
        from rdkit.Chem import Descriptors
        suppl = Chem.SDMolSupplier(str(file_path))
        molecules = [mol for mol in suppl if mol is not None]
        if molecules:
            mol = molecules[0]
            return {
                'Number of molecules': len(molecules),
                'First molecule': mol.GetProp('Name') if mol.HasProp('Name') else 'Unnamed',
                'Molecular weight': f'{Descriptors.MolWt(mol):.2f}',
                'Number of atoms': f'{mol.GetNumAtoms()}',
                'Number of bonds': mol.GetNumBonds(),
                'SMILES': Chem.MolToSmiles(mol)
            }
        else:
            return {
                'Number of molecules': 0,
                'First molecule': 'N/A',
                'Molecular weight': 'N/A',
                'Number of atoms': 'N/A',
                'Number of bonds': 'N/A',
                'SMILES': 'N/A'
            }
    except Exception as e:
        return {
```

```

    'Number of molecules': 'Error',
    'First molecule': 'Error',
    'Molecular weight': 'Error',
    'Number of atoms': 'Error',
    'Number of bonds': 'Error',
    'SMILES': 'Error'
  }

def analyze_xpt_file(file_path):
  """Analyze XPT file and return detailed information."""
  try:
    import pyreadstat
    df, meta = pyreadstat.read_xpt(str(file_path))

    # Get first 5 rows as formatted string
    first_5_rows = df.head(5).to_string(index=True)

    # Get summary statistics
    stats = df.describe().to_string()

    return {
      'Shape': str(df.shape),
      'Columns': str(list(df.columns)),
      'First 5 rows': first_5_rows,
      'Summary statistics': stats
    }
  except Exception as e:
    return {
      'Shape': 'Error: {str(e)}',
      'Columns': 'Error',
      'First 5 rows': 'Error',
      'Summary statistics': 'Error'
    }

def analyze_asnt_file(file_path):
  """Analyze ASNT file and return detailed information."""
  try:
    from Bio import SeqIO
    import xml.etree.ElementTree as ET

    # Check if it's XML
    with open(file_path, 'r') as f:
      first_line = f.readline()

    if first_line.startswith('<xml>'):
      # It's XML, parse as XML
      tree = ET.parse(str(file_path))
      root = tree.getroot()
      records = [root]
    else:
      # Try as ASN.1 sequence
      records = list(SeqIO.parse(str(file_path), "genbank"))

    if records:
      record = records[0]
      if hasattr(record, 'id'): # BioPython record
        return {
          'Number of records': len(records),
          'XML Root': 'N/A (ASN.1 format)',
          'StudyID': 'N/A',
          'Reviewer': 'N/A',
          'AssessmentDate': 'N/A',
          'Findings': 'N/A',
          'Recommendation': 'N/A'
        }
      else: # XML element
        # Extract XML data
        xmldata = {}
        for child in record:
          xmldata[child.tag] = child.text or ''
        return {
          'Number of records': len(records),
          'XML Root': record.tag,
          'StudyID': xmldata.get('StudyID', ''),
          'Reviewer': xmldata.get('Reviewer', ''),
          'AssessmentDate': xmldata.get('AssessmentDate', ''),
          'Findings': xmldata.get('Findings', ''),
          'Recommendation': xmldata.get('Recommendation', '')
        }
    else:
      return {
        'Number of records': 0,
        'XML Root': 'N/A',
        'StudyID': 'N/A',
        'Reviewer': 'N/A',
        'AssessmentDate': 'N/A',
        'Findings': 'N/A',
        'Recommendation': 'N/A'
      }
  except Exception as e:
    return {
      'Number of records': 'Error: {str(e)}',
      'XML Root': 'Error',
      'StudyID': 'Error',
      'Reviewer': 'Error',
      'AssessmentDate': 'Error',
      'Findings': 'Error',
      'Recommendation': 'Error'
    }

# Filter for only the three file types
target_extensions = ['.xpt', '.sdf', '.asnt']

# Create summary data with detailed columns
summary_data = []
for ext in target_extensions:
  if ext in file_extensions:
    files += file_extensions[ext]

for file_path in files:
  file_name = file_path.name

  # Initialize raw data
  row = {
    'File Name': file_name,
    'Extension': ext.upper(),
    'File Path': str(file_path)
  }

  # Add file-specific detailed columns
  if ext == 'sdf':
    sdf_info = analyze_sdf_file(file_path)
    row.update({
      'Number of molecules': sdf_info['Number of molecules'],
      'First molecule': sdf_info['First molecule'],
      'Molecular weight': sdf_info['Molecular weight'],
      'Number of atoms': sdf_info['Number of atoms'],
      'Number of bonds': sdf_info['Number of bonds'],
      'SMILES': sdf_info['SMILES']
    })
  # Add empty columns for other file types
  row.update({
    'Shape': '',
    'Columns': '',
    'First 5 rows': '',
    'Summary statistics': '',
    'Number of records': '',
    'XML Root': '',
    'StudyID': '',
    'Reviewer': '',
    'AssessmentDate': '',
    'Findings': '',
    'Recommendation': ''
  })

  elif ext == '.xpt':
    xpt_info = analyze_xpt_file(file_path)
    row.update({
      'Shape': xpt_info['Shape'],
      'Columns': xpt_info['Columns'],
      'First 5 rows': xpt_info['First 5 rows'],
      'Summary statistics': xpt_info['Summary statistics']
    })
  # Add empty columns for other file types
  row.update({
    'Number of molecules': '',
    'First molecule': '',
    'Molecular weight': '',
    'Number of atoms': '',
    'Number of bonds': '',
    'SMILES': '',
    'Number of records': '',
    'XML Root': '',
    'StudyID': '',
    'Reviewer': '',
    'AssessmentDate': '',
    'Findings': '',
    'Recommendation': ''
  })

  elif ext == '.asnt':
    asnt_info = analyze_asnt_file(file_path)
    row.update({
      'Number of records': asnt_info['Number of records'],
      'XML Root': asnt_info['XML Root'],
      'StudyID': asnt_info['StudyID'],
      'Reviewer': asnt_info['Reviewer'],
      'AssessmentDate': asnt_info['AssessmentDate'],
      'Findings': asnt_info['Findings'],
      'Recommendation': asnt_info['Recommendation']
    })
  # Add empty columns for other file types
  row.update({
    'Number of molecules': '',
    'First molecule': '',
    'Molecular weight': '',
    'Number of atoms': '',
    'Number of bonds': '',
    'SMILES': '',
    'Shape': '',
    'Columns': '',
    'First 5 rows': '',
    'Summary statistics': ''
  })

  summary_data.append(row)

```

# Create DataFrame and display																			
summary_df = pd.DataFrame(summary_data)																			
summary_df																			
[16:12:83] ER80R: Atom line too short: ' -0.0015 1.2895 8.8900 C' on line 5																			
[16:12:83] ER80R: moving to the beginning of the next molecule																			
Out [7]:																			
File Name	Extension	File Path	Shape	Columns	First 5 rows	Summary statistics	Number of molecules	First molecule	Molecular weight	Number of atoms	Number of bonds	SMILES	Number of records	XML Root	StudyID	Reviewer	AssessmentDate	Findings	Recommendation
0	lb.xpt	.XPT	reports/study1234/datasets/dat...	(2, 5)	["STUDYID", "USUBJID", "LBTEST", "LBSTRESN", "LBTR..."]	STUDYID: 2.000000 mean 51.10...	LBSTRESN count 51.10...												
1	ae.xpt	.XPT	reports/study1234/datasets/dat...	(2, 5)	["STUDYID", "USUBJID", "AETERM", "AESEV", "AER..."]	STUDYID: 2.000000 mean 48.00...	USUBJID count 48.00...	AETERM: AEEV count 48.00...	AEEV: AER...										
2	dm.xpt	.XPT	reports/study1234/datasets/dat...	(2, 4)	["STUDYID", "USUBJID", "AGE", "SEX"]	STUDYID: 2.000000 mean 48.00...	USUBJID: AGE count 48.00...	AGE count 48.00 mean 48.00...	SEX p0 ABC123	SUBJ0...									
3	compound.sdf	.SDF	reports/study1234/datasets/dat...				0	N/A	N/A	N/A	N/A	N/A							
4	Structure2D_COMPOUND_CID_197365.sdf	.SDF	reports/study1234/datasets/dat...				1	197365	699.99	44	46	c4ccc5nc(CCCc6cc[N]CCO)C(Cl)...	CN1CCN(c2ccc3nc ...						
5	Conformer3D_COMPOUND_CID_197366.sdf	.SDF	reports/study1234/datasets/dat...				1	197366	590.60	41	46	c4ccc5nc(CCCc6cc[N]CCO)C(Cl)...	CN1CCN(c2ccc3nc ...						
6	Conformer3D_COMPOUND_CID_197366.asnt	ASNT	reports/study1234/datasets/dat...										0	N/A	N/A	N/A	N/A	N/A	
7	Structure2D_COMPOUND_CID_197365.asnt	ASNT	reports/study1234/datasets/dat...										0	N/A	N/A	N/A	N/A	N/A	
8	assessment.asnt	ASNT	reports/study1234/datasets/dat...										1	AssessmentTemplate	STUDY1234	Dr. Smith	2025-11-05	\\n Approved for next review phase	

Save Summary Table to CSV

```
In [8]: # Save the summary table to CSV
csv_filename = 'SD_Study-Data-files.csv'
summary_df.to_csv(csv_filename, index=False)
print(f"Summary table saved to {csv_filename}")
print(f"CSV file contains {len(summary_df)} rows and {len(summary_df.columns)} columns")
print(f"Columns: {list(summary_df.columns)}")
Summary table saved to SD_Study-Data-files.csv
CSV file contains 9 rows and 28 columns
Columns: ['File Name', 'Extension', 'File Path', 'Shape', 'Columns', 'First 5 rows', 'Summary statistics', 'Number of molecules', 'First molecule', 'Molecular weight', 'Number of atoms', 'Number of bonds', 'SMILES', 'Number of records', 'XML Root', 'StudyID', 'Reviewer', 'AssessmentDate', 'Findings', 'Recommendation']
```